

AtriCure®**Isolator™ Transpolar™ pen****Instructions for Use**

The Isolator™ Transpolar™ pen System is comprised of the AtriCure® Ablation and Sensing Unit (ASU), the Isolator™ Transpolar™ pen (Pen) or Isolator™ Long Pen TT (Pen), Footswitch, and ASU Source Switch. The Pen is a single patient use electrocautery instrument designed for use only with the ASU/ASU2. The Pen is used to ablate soft tissues. When the Pen is connected to the ASU, the ASU provides the bipolar radiofrequency (RF) energy flowing between both electrodes of the Pen. The Operator controls the application of this RF energy by pressing the Footswitch.

ISOLATOR™ TRANSPOLAR™ PEN ILLUSTRATION AND NOMENCLATURE

1. Distal Tip
2. Electrodes
3. Malleable Shaft
4. Handle

ISOLATOR™ TT LONG PEN ILLUSTRATION AND NOMENCLATURE

1. Distal Tip
2. Electrodes
3. Malleable Shaft
4. Rigid Shaft

INDICATION FOR USE

The Isolator™ Transpolar™ pen is a sterile, single use electrocautery device intended for the ablation of soft tissues using bipolar radiofrequency (RF) energy during General, ENT, Thoracic, Gynecology, and Urology surgical procedures.

CONTRAINDICATIONS

- The device is not intended for contraceptive tubal coagulation (permanent female sterilization).
- Do not ablate directly on cardiac valves.

POTENTIAL COMPLICATIONS

Possible complications related to the creation of spot or linear lesions in cardiac and soft tissues are:

- Tissue perforation
- Postoperative embolic complications
- Extension of extracorporeal bypass
- Perioperative heart rhythm disturbance (atrial and/or ventricular)
- Pericardial effusion and tamponade
- Injury to the great vessels
- Valve leaflet damage
- Conduction disturbances (SA/AV node)
- Acute ischemic myocardial event

WARNINGS

Do not touch the electrodes of the Pen while activating the ASU. Touching the Pen electrodes during ASU activation could result in an electrical shock or burn to the operator.

- Do not touch the electrodes of the Pen or metal staples or clips, or to sutures while activating the ASU. This may damage the Pen or tissue, or result in an incomplete ablation.
- Do not use abrasive cleaners or electrocautery tip cleaners to clean debris from the Distal Tip. Do not turn on ASU while cleaning the device. Use of abrasive cleaners or electrocautery tip cleaners can damage the electrodes and result in device failure. Use saline-soaked gauze for cleaning debris.
- Do not immerse the Pen in liquids as this may damage the device.
- To avoid shock/burn hazard, always wear the appropriate surgical gloves when using the Pen and ASU.
- To avoid shock/burn hazards, always remove the Pen from the patient during defibrillation.
- No representation or warranty is made that failure or cessation of function of the device will not result in an adverse event or that medical complication (including tissue perforation) will not follow the procedure, or that the use of the device will in all cases restore adequate cardiac function.


PRECAUTIONS

- Read all instructions carefully for the AtriCure® ASU Isolator™ Transpolar™ pen, ASU Source Switch, and any auxiliary device being used prior to using the device. Failure to properly follow instructions may lead to electrical or thermal injury and may result in improper functioning of the device.
- Use of the Pen should be limited to properly trained and qualified medical personnel. Proper surgical procedures and techniques are the responsibility of the medical professional. Each surgeon must evaluate the appropriateness of any procedure based on their own medical training and experience, and the type of surgical procedure.
- Patient and procedure selection is solely a medical responsibility and the outcome is dependent on many variables including patient pathology and surgical and perfusion procedures.
- Variations in specific procedures may occur due to individual physician techniques and patient anatomy.
- To avoid damage to the device, do not drop or toss the Pen. If the Pen is dropped, do not use. Replace with a new Pen.
- Do not use the Pen in the presence of flammable materials.
- Do not re-sterilize or reuse the Pen. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another.
- The distal tip of the Pen must be kept clean of debris during surgery to avoid loss of power. Before activating the ASU, inspect the area at the distal tip of the Pen for foreign matter. Foreign matter captured on the tip will adversely affect the ablation.
- The Pen is only compatible with the AtriCure® ASU2 and ASU3. Use of the Pen with another manufacturer's generator may damage the device and result in patient injury.
- The ASU should only be used with approved cardiac pacing and sensing devices.
- Do not use the Pen for coagulation or ablation of veins or arteries.
- Before activating the ASU, inspect the area at the distal tip of the Pen for foreign matter. Foreign matter captured on the tip will adversely affect the ablation.
- Use caution to avoid trauma to tissues not within the target area of ablation. Tissue and/or structures behind the targeted tissue should be protected from potential thermal spread.
- The Pen has an eight hour useful life that is tracked by the ASU. If the Pen no longer functions and the ASU will display a message indicating that the Pen must be replaced.
- It is the responsibility of the user to dispose of this device in accordance with local regulations.
- Excessive bending of the malleable stainless steel shaft will cause the shaft to harden and may increase the potential for breakage.

INSTRUCTIONS FOR USE

1. Using sterile technique, remove the Pen from its packaging.

2. With the Connector Alignment Arrow symbol in the 12 o'clock position, push the Connector into the Pen receptacle on the front of the ASU. Verify that the connections between the Pen and the ASU are secure. If the connections are loose, do not use the Pen. Inspect the Cable and do not use the Pen if the cable is frayed or the insulation is damaged.

**3. Painting ablation technique**

- Maintaining visualization, move the distal tip gently across the targeted cardiac tissue.
- While maintaining continuous contact between the tissue and the electrodes, move the device continuously in an oscillating manner at a rate of approximately 1cm/sec.

| |
|----------------------------|
| Painting Lesion Thickness* |
| 20 seconds |
| 2 cm oscillation |
| 2.0 - 4.0 mm |

*Data was obtained from ablations performed on excised bovine myocardium and represent 95% confidence intervals. Results may vary based on tissue properties.

4. Stamping Ablation technique:

- Apply constant firm pressure to the tissue without movement. Maintain full contact of the electrode surface with the tissue. A stamping lesion is approximately 8 mm x 6 mm.
- If creating longer lesions with the stamping technique, overlap the contiguous ablations by 50% to ensure a continuous and complete lesion.

| |
|-------------------------|
| Stamp Lesion Thickness* |
| 10 seconds |
| 15 seconds |
| 3.3 - 3.8 mm |
| 3.8 - 4.4 mm |

*Data was obtained from ablations performed on excised bovine myocardium and represent 95% confidence intervals. Results may vary based on tissue properties.

5. Press the Footswitch to activate the ASU.

6. When the Footswitch is pressed, the ASU will emit an audible tone indicating that current is flowing between the electrodes located at the distal tip of the Pen and through the tissue.

7. Inspect the surgical area to ensure adequate ablation.

Between ablations, wipe the distal tip clean with a saline-soaked gauze pad. Important: For optimal performance, keep the Pen electrodes clear of coagulum. To ensure the electrodes are clear of coagulum:

- Use saline-soaked gauze to clean the electrodes after each ablation. The coagulum is much easier to remove within the first several seconds after ablation. In a brief period of time, the coagulum may dry making removal of coagulum more difficult.
- Check both electrodes before each ablation to ensure that the gold of each electrode is visible and coagulum is removed.

9. If the Pen is idle between ablations, place the Pen distal tip onto saline-soaked gauze to prevent any coagulum not cleaned off the electrodes from drying. Repeat ablation if necessary.

10. Upon completion of the surgical procedure, disconnect the Pen from the generator and discard the Pen after use. Follow local governing ordinances and recycling plans regarding disposal or recycling of device component.

HOW SUPPLIED

The Isolator™ Transpolar™ pen is supplied as a STERILE instrument and is for single patient use only. Sterility is guaranteed unless the package is opened or damaged. Do not resterilize.

RETURN OF USED PRODUCT

If any return of this product must be returned to AtriCure®, Inc., a return goods authorization (RGA) number is required from AtriCure®, Inc., prior to shipping. If the product has been in contact with body or blood fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton, to prevent damage during shipment; and it should be properly labeled with an RGA number and an indication of the biohazardous nature of the contents of shipment. Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number may be obtained from AtriCure®, Inc.

ATTENTION! It is the responsibility of the health care institution to adequately prepare and identify the products for shipment.

DISCLAIMER STATEMENTS

Users assume responsibility for approving the acceptable condition of this product before it is used, and for ensuring that the product is only used in the manner described in these instructions for use, including, but not limited to, ensuring that the product is not re-used.

Under no circumstances will AtriCure, Inc. be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse or re-use of this product, including any loss, damage, or expense which is related to personal injury or damage to property.

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