

**US Indications:** AtriCure's cryoICE® cryoSPHERE™ cryoablation probes are sterile, single use devices intended for use in blocking pain by temporarily ablating peripheral nerves.

Please review the Instructions for Use for a complete listing of contraindications, warnings, precautions and potential adverse events prior to using these devices.

**Rx Only.**

## *cryo*SPHERE™

- Enhanced Flexibility
- Uniform Application
- Engineered for Intercostal Nerve Cryoablation

### **Indications**

AtriCure's cryoICE® cryoSPHERE™ cryoablation probes are sterile, single use devices intended for use in blocking pain by temporarily ablating peripheral nerves.

### **ATRICURE, INC.**

7555 Innovation Way  
Mason, Ohio 45040 USA  
Customer Service: +1 (866) 349-2342  
CustomerService@AtriCure.com  
[www.AtriCure.com](http://www.AtriCure.com)

PM-US-0058B-0121-G

# *cryo*SPHERE™


The Coolest Innovation  
in Pain Management.



**AtriCure**

**AtriCure**


# The Coolest Innovation in Pain Management.



## Features

<b>8mm Ball Tip</b>	Designed for minimally invasive applications and placement on intercostal nerve
<b>Bendable Distal Shaft</b>	Bendable distal shaft up to 180 degrees
<b>Braided Insulation</b>	Technology that is designed to maintain a temperature warmer than -10°C in tissue
<b>Reinforced Proximal Shaft</b>	Provides rigidity at the level of the incision points
<b>Ergonomic Handle Design</b>	Provides ease of use and comfort during cryoablation
<b>Active Defrost</b>	Allows probe to be easily removed from cryoablation site

## Malleable Cryoablation Probe



## cryoICE BOX® - V6

ACM1 (US)



## cryoICE Cart & Upgrade Kits

ASCS (cart with double tank holder, US only)

ASCC (single tank holder upgrade)

ASCC2 (double tank holder upgrade)

20 lbs N<sub>2</sub>O tanks not provided by AtriCure



## Foot Switch (optional)

For use with ACM1/ACM2  
FTSWITCH

